



12098

Black back rest and seat cushion for Extra chair



FEATURES

Mattress for sedan chair designed to increase patients' comfort and to protect the seat frame at the same time. Made of expanded polyurethane at slow memory, it has a fire-proof and anti-tear cover that allows it to be easily removable and washed. Designed for Extra sedan chair, it is supplied with Velcro® fixing strips that make it suitable for the entire range of sedan chairs. Available in orange and black colours.

Dimensions:

SEAT 42 x 43,5 x 6 [cm] (HxLxD)

BACK REST 46,5 x 43,5 x 3 [cm] (HxLxD)

*All Pictures shown are for illustration purpose only.

ACCESSORY OF



Extra 670/BR
Track assisted stair chair with arms



678 Extra Lite Basic
Extra Lite Basic, foldable carrying chair



Extra Lite Standard 676
Foldable carrying chair with 4 wheels, footrest and arms



Extra Ergolift 671
Powered stair chair