



7270 PROOF MBS-01 MEBER BARIATRIC STRETCHER 7270

Self loading stretcher for bariatric transport



FEATURES

Self loading stretcher for bariatric transport made of a special high resistance stainless steel alloy. The frame is reinforced to ensure a 350 kg maximum payload compared to others MeBer self loading stretchers. It is supplied with 2 frontal Ø125x25 and 2 back Ø125x25 wheels for maximum safety during loading operations. The increased 4 Ø200x50 mm wheels, 2 fixed with brake and 2 rotating with self-centring system, allow easy manoeuvrability even with heavy loads. The rigid platform, composed by reinforced thermoplastic material, is easy to clean and disinfect. It also has an adjustable backrest with lockable gas spring for inclinations from 0° to 80° and fully folding sides. Vertical legs are protected by plastic buffers. This stretcher is supplied complete with mattress (art. 7013) and belts (N. 2 art. 608-A/MEB or 608-N/MEB and n. 1 art. 612/AN-2/MEB or art. 612/N-2/MEB).

Thanks to the increased load capacity it can be used in coupling with the MBP bariatric platform (Article code 12080) and bariatric wings set (art. 12084).



MAIN FEATURES

Usage	Transport
Lying surface	Without Trendelenburg
Type	Self loading stretcher
Height (cm)	92
Width (cm)	56
Weight (kg)	47,5
Loading capacity (kg)	350
Folded height (cm)	25
Max loading height (cm)	64,5
Length (cm)	197
Material	Stainless steel
Bariatric	

IMAGES



ACCESSORIES



12025

Black multi usage object holder for automatic loading stretcher



12084

Bariatric wings set for self-loading stretchers complete with mattresses and belts with reel



608/A-MEB

2 pieces orange belt with metallic hook and MeBer band (L. 155 cm)



608/N-MEB

2 pieces black belt with metallic hook and MeBer band (L. 155 cm)



612/AN-2/MEB

Adjustable thoracic abdominal belt orange/black with MeBer band and 4 points fastening for self loading stretchers and sedan chairs



612/N-2/MEB

Adjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self loading stretchers and sedan chairs



7013/N

3 pcs. Skay fire-proof electrowelded mattress



12106

Mattress for AIR-COM stretchers - Anti-bedsore mattress with alternate segments



12070/BG

"COVER OX 2-3 BG"- professional cylinder holder for 2 and 3 lt oxygen bottles



12070/BR

COVER OX2-3 BR professional cylinder holder with belt and orange MeBer logo for 2 and 3 lt oxygen bottles



12072/BG

"COVER OX 5 BG"- professional cylinder holder for 5 lt oxygen bottles



12072/BR

"COVER OX 5 BR"- professional cylinder holder for 5 lt oxygen bottles



938

Oximeb oxygen cylinder holder



12054

Oximeb Maxi - Oxygen cylinder holder



12154

"BEAR PODS" Bariatric containment system for self-loading stretchers



707

Red 6-handle fireproof cover with Me.Ber tape for self-loading stretcher



955 PROOF

EN 1789 fixing system for Proof roll-in stretchers



12080-001

Monting kit for bariatric platform on self loading stretcher



12080

"MPB" wide bariatric platform for



7015

Folding sky mattress black for bariatric MBP platform

self-loading stretchers Frog
Mercury and Gecko complete of
belts with reel and fixing for
stretcher

STANDARD ACCESSORIES



7013/N

3 pcs. Skay fire-proof electrowelded
mattress



608/A-MEB

2 pieces orange belt with metallic
hook and MeBer band (L. 155 cm)



612/AN-2/MEB

Adjustable thoracic abdominal belt
orange/black with MeBer band and
4 points fastening for self loading
stretchers and sedan chairs