



## 7270 PROOF MBS-01 MEBER BARIATRIC STRETCHER 7270

---

Self loading stretcher for bariatric transport



---

### FEATURES

Self loading stretcher for bariatric transport made of a special high resistance stainless steel alloy. The frame is reinforced to ensure a 350 kg maximum payload compared to others MeBer self loading stretchers. It is supplied with 2 frontal Ø125x25 and 2 back Ø125x25 wheels for maximum safety during loading operations. The increased 4 Ø200 mm wheels, 2 fixed with brake and 2 rotating with self-centring system, allow easy manoeuvrability even with heavy loads. The rigid platform, composed by aluminium, is easy to clean and disinfect. It also has an adjustable backrest with lockable gas spring for inclinations from 0° to 80° and fully folding sides. Vertical legs are protected by plastic buffers. This stretcher is supplied complete with mattress (art. 7013) and belts (N. 2 art. 608-A/MEB or 608-N/MEB and n. 1 art. 612/AN-2/MEB or art. 612/N-2/MEB).

Thanks to the increased load capacity it can be used in coupling with the MBP bariatric platform (Article code 12080) and bariatric wings set (art. 12084).



---

## MAIN FEATURES

---

Usage	Transport
Lying surface	Without Trendelenburg
Type	Self loading stretcher
Height (cm)	92
Width (cm)	56
Weight (kg)	49,5
Loading capacity (kg)	350
Folded height (cm)	25
Max loading height (cm)	64,5
Length (cm)	197
Material	Stainless steel
Bariatric	

---

## IMAGES

---



## ACCESSORIES



**12025**

Black multi usage object holder for automatic loading stretcher



**12084**

Bariatric wings set for self-loading stretchers complete with mattresses and belts with reel



**608/A-MEB**

2 pieces orange belt with metallic hook and MeBer band (L. 155 cm)



**608/N-MEB**

2 pieces black belt with metallic hook and MeBer band (L. 155 cm)



**612/AN-2/MEB**

Adjustable thoracic abdominal belt orange/black with MeBer band and 4 points fastening for self loading stretchers and sedan chairs



**612/N-2/MEB**

Adjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self loading stretchers and sedan chairs



**7013/N**

3 pcs. Skay fire-proof electrowelded mattress



**12106**

Mattress for AIR-COM stretchers - Anti-bedsore mattress with alternate segments



**12070/BG**

"COVER OX 2-3 BG"- professional cylinder holder for 2 and 3 lt oxygen bottles



**12070/BR**

COVER OX2-3 BR professional cylinder holder with belt and orange MeBer logo for 2 and 3 lt oxygen bottles



**12072/BG**

"COVER OX 5 BG"- professional cylinder holder for 5 lt oxygen bottles



**12072/BR**

"COVER OX 5 BR"- professional cylinder holder for 5 lt oxygen bottles



**938**

Oximeb oxygen cylinder holder



**12054**

Oximeb Maxi - Oxygen cylinder holder



**12154**

"BEAR PODS" Bariatric containment system for self-loading stretchers



**707**

Red 6-handle fireproof cover with Me.Ber tape for self-loading stretcher



**955 PROOF**

EN 1789 fixing system for Proof roll-in stretchers



**12080-001**

Monting kit for bariatric platform on self loading stretcher



**12080**

"MPB" wide bariatric platform for



**7015**

Folding sky mattress black for bariatric MBP platform

self-loading stretchers Frog  
Mercury and Gecko complete of  
belts with reel and fixing for  
stretcher

---

---

## STANDARD ACCESSORIES

---



**7013/N**

3 pcs. Skay fire-proof electrowelded  
mattress



**608/A-MEB**

2 pieces orange belt with metallic  
hook and MeBer band (L. 155 cm)



**612/AN-2/MEB**

Adjustable thoracic abdominal belt  
orange/black with MeBer band and  
4 points fastening for self loading  
stretchers and sedan chairs

---