



1398

White/blue absorbing and transfer sheet cm. 220x100



FEATURES

Transfer sheet with a 3-layer construction (PE-foil, nonwoven and cellulose fluff), it is highly resistant (up to 150 kg) and it allows the absorbency of fluids and exudates (up to 3,5 liters).

It is an effective and hygienic solution, developed with the goal of favoring the transfer of patients, where other methods of transfer could be difficult, or for being used in complement to other means of recovery and transport (vacuum mattresses, rescue stretchers, scoop stretchers, Toboga stretchers etc.); after its use the sheet can be easily eliminated and the carriers or stretchers are ready for the next use.

The innovative sheet 2 in 1 is a disposable device of class I (according to the European Directive 93/42/CEE).

Sheet dimensions: 220x100 cm

Dimensions of the absorbing part: 200x50 cm

Weight: 350g

ATTENTION: The application of the transfer sheet should be always according to the application requirements/indications on its package.

MAIN FEATURES

Width (cm)	100
Length (cm)	220
Loading capacity (kg)	150

IMAGES



SEE ALSO



15180
Big Gurney
