



12106

Mattress for AIR-COM stretchers - Anti-bedsores mattress with alternate segments



FEATURES

Upper mattress made in PVC and internally composed by 18 independent, interchangeable air cells. Particularly suited for the positioning of patients with a risk of decubitus ulcers. AIR-COM can be positioned on top of any type of stretchers mattress. The underlying part, made of anti-slip textile, is completed by 3 belt loops to fix it to the stretcher.

Dimensions: cm 180x46x10

*All Pictures shown are for illustration purpose only

MAIN FEATURES

Length (cm)	180
Width (cm)	46
Height (cm)	10

IMAGES



ACCESSORIES



Art. 12108
Compressor for AIR-COM mattress

ACCESSORY OF



Mercury Cinque 7077/4RG PROOF
Certified 5 levels self-loading
stretcher with 4 swivel wheels -
high version



**Mercury Lite Cinque 7097/4RG
Proof**
Certified aluminium 5 levels self-
loading stretcher with 4 swivel
wheels - high version



Mercury Lite 7082/4RG Proof
Certified aluminium variable height
stretcher high version with 4
rotating wheels



Mercury Cinque 7075/4RG PROOF
Certified 5 levels self-loading
stretcher with 4 swivel wheels



**Mercury Lite Cinque 7095/4RG
Proof**
Certified aluminium 5 levels self-
loading stretcher with 4 swivel
wheels



Mercury EL 7092 PROOF
Certified self-loading stretcher



Mercury EL 7094 PROOF
Certified self-loading stretcher -
high version



7265 Proof
Aluminium self loading stretcher
"Frog lite Advance" with
extractable handles and
Trendelenburg
