



## 12106

Mattress for AIR-COM stretchers - Anti-bedsore mattress with alternate segments



### FEATURES

Upper mattress made in PVC and internally composed by 18 independent, interchangeable air cells. Particularly suited for the positioning of patients with a risk of decubitus ulcers. AIR-COM can be positioned on top of any type of stretchers mattress. The underlying part, made of anti-slip textile, is completed by 3 belt loops to fix it to the stretcher.

**Dimensions: cm 180x46x10**

\*All Pictures shown are for illustration purpose only

### MAIN FEATURES

Length (cm)	180
Width (cm)	46
Height (cm)	10

### IMAGES



---

ACCESSORIES

---



**Art. 12108**  
Compressor for AIR-COM mattress

---

ACCESSORY OF

---



**Mercury Cinque 7077/4RG PROOF**  
Certified 5 levels self-loading  
stretcher with 4 swivel wheels -  
high version



**Mercury Lite Cinque 7097/4RG  
Proof**  
Certified aluminium 5 levels self-  
loading stretcher with 4 swivel  
wheels - high version



**Mercury Lite 7082/4RG Proof**  
Certified aluminium variable height  
stretcher high version with 4  
rotating wheels



**Mercury Cinque 7075/4RG PROOF**  
Certified 5 levels self-loading  
stretcher with 4 swivel wheels



**Mercury Lite Cinque 7095/4RG  
Proof**  
Certified aluminium 5 levels self-  
loading stretcher with 4 swivel  
wheels



**Mercury EL 7092 PROOF**  
Certified self-loading stretcher



**Mercury EL 7094 PROOF**  
Certified self-loading stretcher -  
high version



**7265 Proof**  
Aluminium self loading stretcher  
"Frog lite Advance" with  
exctractable handles and  
Trendelenburg