



7034


Fixed trolley for incubators



FEATURES

Fixed height incubator trolley made of stainless steel structure with No. 4 wheels of \varnothing 200 mm (2 with brake) for transport inside hospitals. Compact dimensions for easy manoeuvrability. This trolley is equipped with a universal interface for anchoring the incubator.


MAIN FEATURES

Pediatric	
Material	Stainless steel
Height (cm)	74
Width (cm)	57
Length (cm)	90
Weight (kg)	30,5

SEE ALSO



12050
Extended trolley incubator
interface



12048
Certified interface for incubator
fixation for self loading trolleys

