



## 7240/INC PROOF Frog

"FROG PLUS INC PROOF" self loading trolley for incubators  
predisposed for certified fixing system



### FEATURES

Made with predisposition for our certified fixing device Art. 955 Proof, it is perfect for transport of incubators on rescue vehicles. Designed without back rest, side protections and rigid platform. This stretcher is part of "FROG PLUS" self loading stretchers range.



Device certified according to the European harmonized safety standards UNI EN1865 - UNI EN 1789.

### MAIN FEATURES

Pediatric	
Height (cm)	82
Max loading height (cm)	66
Min loading height (cm)	60
Folded height (cm)	27
Width (cm)	56
Length (cm)	195
Length with folded trolley (cm)	166
Weight (kg)	38,4
Loading capacity (kg)	250

---

## SEE ALSO

---



**Frog 7240 PROOF**  
"FROG PLUS" self loading stretcher  
with trendelenburg certified

---

---

## ACCESSORIES

---



**12050**  
Extended trolley incubator  
interface



**12048**  
Certified interface for incubator  
fixation for self loading trolleys



**955 PROOF**  
EN 1789 fixing system for Proof roll-  
in stretchers

---