



1510 PROOF Fixed stretcher

Fixed stretcher with nylon cover



FEATURES

Fixed stretcher built in anodized aluminium, canvas in fire-resistant and easily cleanable nylon. This model is suitable for our ambulance trolleys (Titanus) or as a hand-carried stretcher.

- 4 extractable and adjustable handles
- 4 wheels diam. 100mm
- Adjustable backrest
- Supplied with straps (art. 604/MEB)
- Foldable side rails

Supplied with belts only in combination with 7060/4RG trolley:

- n. 2 pcs. art. 608-A/MEB o 608-N/MEB
- n. 1 pc art. 612/AN-2/MEB o art. 612/N-2/MEB



Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789

MAIN FEATURES

Size	Adults
Usage	Transport
Lying surface	Without Trendelenburg
Height (cm)	16
Width (cm)	57
Length (cm)	196
Max length (cm)	233
Weight (kg)	15
Loading capacity (kg)	170
Wheel type	4 Small Fixed Wheels
10G Certification	
Material	Anodized aluminium Nylon

IMAGES



SEE ALSO



Fixed stretcher 1515 PROOF
Fixed stretcher with rigid bed



Titanus 7060/4RG PROOF
Certified stainless steel variable
height trolley with four rotating
wheels



Titanus 7062/4RG PROOF
Stainless steel variable height
trolley with four rotating wheels
high version

ACCESSORIES



936/A

Phlebotomy stand for stretchers
Art. 1510 Proof and Art. 1515 Proof



608/A-MEB

2 pieces orange belt with metallic
hook and MeBer band (L. 155 cm)



608/N-MEB

2 pieces black belt with metallic
hook and MeBer band (L. 155 cm)



612/AN-2/MEB

Adjustable thoracic abdominal belt
orange/black with MeBer band and
4 points fastening for self loading
stretchers and sedan chairs



612/N-2/MEB

Adjustable thoracic abdominal belt
black with MeBer band with 4
points fastening for self loading
stretchers and sedan chairs

STANDARD ACCESSORIES



604/MEB

Orange bag with 2 straps with
plastic buckle and ME.BER. ribbon