



## 1510 PROOF Fixed stretcher

---

Fixed stretcher with nylon cover



---

### FEATURES

Fixed stretcher built in anodized aluminium, canvas in fire-resistant and easily cleanable nylon. This model is suitable for our ambulance trolleys (Titanus) or as a hand-carried stretcher.

- 4 extractable and adjustable handles
- 4 wheels diam. 100mm
- Adjustable backrest
- Supplied with straps (art. 604/MEB)
- Foldable side rails

Supplied with belts only in combination with 7060/4RG trolley:

- n. 2 pcs. art. 608-A/MEB o 608-N/MEB
- n. 1 pc art. 612/AN-2/MEB o art. 612/N-2/MEB



Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789

---

## MAIN FEATURES

---

Size	Adults
Usage	Transport
Lying surface	Without Trendelenburg
Height (cm)	16
Width (cm)	57
Length (cm)	196
Max length (cm)	233
Weight (kg)	15
Loading capacity (kg)	170
Wheel type	4 Small Fixed Wheels
10G Certification	
Material	Anodized aluminium Nylon

---

## IMAGES

---



---

## SEE ALSO

---



**Fixed stretcher 1515 PROOF**  
Fixed stretcher with rigid bed



**Titanus 7060/4RG PROOF**  
Certified stainless steel variable  
height trolley with four rotating  
wheels



**Titanus 7062/4RG PROOF**  
Stainless steel variable height  
trolley with four rotating wheels  
high version

---

---

## ACCESSORIES

---



**936/A**

Phlebotomy stand for stretchers  
Art. 1510 Proof and Art. 1515 Proof



**608/A-MEB**

2 pieces orange belt with metallic  
hook and MeBer band (L. 155 cm)



**608/N-MEB**

2 pieces black belt with metallic  
hook and MeBer band (L. 155 cm)



**612/AN-2/MEB**

Adjustable thoracic abdominal belt  
orange/black with MeBer band and  
4 points fastening for self loading  
stretchers and sedan chairs



**612/N-2/MEB**

Adjustable thoracic abdominal belt  
black with MeBer band with 4  
points fastening for self loading  
stretchers and sedan chairs

---

## STANDARD ACCESSORIES

---



**604/MEB**

Orange bag with 2 straps with  
plastic buckle and ME.BER. ribbon