



980

"MEBER-X" Multi level cot fastening system



FEATURES

"MEBER-X" Multi level cot fastening system composed of front and side fixations and rear safety hook.

Front fixation

Length: 540 mm

Width: 680 mm

Height: 140 mm

Side fixation

Length: 750 mm

Width: 80 mm

Height: 240 mm

Rear safety hook

Length: 100 mm

Width: 20 mm

Height: 50 mm

ACCESSORY OF



MeBer-X 7300
"MeBer-X" Multi-level stretcher
complete with mattress and belts
