



12050

Extended trolley incubator interface



FEATURES

Extended trolley incubator interface

MAIN FEATURES

Length (cm)	135
Width (cm)	57
Height (cm)	2,5
Weight (kg)	20
Pediatric	

SEE ALSO



12048
Certified interface for incubator
fixation for self loading trolleys



Self loading trolleys for incubator



7034
Fixed trolley for incubators

ACCESSORY OF



Frog 7240/INC PROOF
"FROG PLUS INC PROOF" self
loading trolley for incubators
predisposed for certified fixing
system



Frog Lite 7260/INC PROOF
Self loading stretcher for incubators
predisposed for certified fixing
system



Winner 910/GIN PROOF
Self loading trolley for incubators
predisposed for certified fixing
system



Winner 910/GIN
Self loading trolley for incubators
