



## 15146

---

Anatomical anti-bedsore cushion in polyurethane foam with polyurethane gel insert



---

### FEATURES

Anatomically shaped anti-bedsore cushion, composed of a "slow memory" polyurethane foam base. The device is covered with a removable cover which can be machine washed at 30°. The upper part of the cover in contact with the patient is elastic and hypoallergenic, the lower part is anti-slip. Ideal for medium risk decubitus.

Technical specifications:

Material: slow memory polyurethane foam (density 52-58 kg/m<sup>3</sup>)

Cover material: single spread polyurethane in the upper part and PVC polyester in the lower part

Dimensions: cm 41x41x7 H

Characteristics: removable cover


---

### SEE ALSO



**15148**

Anti-bedsore cushion in polyurethane gel with expanded foam base



**15142**

Anti-bedsore fluid gel cushion



**15144**

Viscoelastic anti-bedsore cushion



**15150**

Anti-bedsore air cushion with connected cells