



## 12018N

---

Telescopic foldable phleboclysis stand for self loading stretchers



---

### FEATURES

"P-120" Aluminium telescopic foldable phleboclysis stand for self loading stretchers.

Dimensions: min. height cm 45 - max. height cm 128

Height 1st stage: cm 70

Height 2nd stage: cm 98

---

### IMAGES



---

## SEE ALSO

---



**936**

Phlebotomy stand for automatic loading stretcher Winner and Frog



**936/A**

Phlebotomy stand for stretchers Art. 1510 Proof and Art. 1515 Proof

---

## ACCESSORY OF

---



**Mercury Lite 7080/4RG Proof**

Ambulance cot in aluminium with variable height and four swivel castors



**Mercury 7070/4RG PROOF**

Stainless steel variable height stretcher certified with 4 rotating wheels



**Mercury 7072/4RG PROOF**

Certified stainless steel variable height stretcher high version with 4 rotating wheels



**Mercury Cinque 7075/4RG PROOF**

Certified 5 levels self-loading stretcher with 4 swivel wheels



**Mercury Cinque 7077/4RG PROOF**

Certified 5 levels self-loading stretcher with 4 swivel wheels - high version



**Mercury Lite Cinque 7095/4RG Proof**

Certified aluminium 5 levels self-loading stretcher with 4 swivel wheels



**Mercury Lite Cinque 7097/4RG Proof**

Certified aluminium 5 levels self-loading stretcher with 4 swivel wheels - high version



**Mercury EL 7092 PROOF**  
Certified self-loading stretcher



**Mercury EL 7094 PROOF**

Certified self-loading stretcher - high version



**MeBer-X 7300**

"MeBer-X" Multi-level stretcher complete with mattress and belts