



12018N

Telescopic foldable phleboclysis stand for self loading stretchers



FEATURES

"P-120" Aluminium telescopic foldable phleboclysis stand for self loading stretchers.

Dimensions: min. height cm 45 - max. height cm 128

Height 1st stage: cm 70

Height 2nd stage: cm 98

IMAGES



SEE ALSO



936

Phlebotomy stand for automatic loading stretcher Winner and Frog



936/A

Phlebotomy stand for stretchers Art. 1510 Proof and Art. 1515 Proof

ACCESSORY OF



Mercury Lite 7080/4RG Proof

Ambulance cot in aluminium with variable height and four swivel castors



MeBer-X 7300

"MeBer-X" Multi-level stretcher complete with mattress and belts



Mercury Cinque 7075/4RG PROOF

Certified 5 levels self-loading stretcher with 4 swivel wheels



Mercury Cinque 7077/4RG PROOF

Certified 5 levels self-loading stretcher with 4 swivel wheels - high version



Mercury Lite Cinque 7095/4RG Proof

Certified aluminium 5 levels self-loading stretcher with 4 swivel wheels



Mercury Lite Cinque 7097/4RG Proof

Certified aluminium 5 levels self-loading stretcher with 4 swivel wheels - high version



Mercury EL 7092 PROOF

Certified self-loading stretcher



Mercury EL 7094 PROOF

Certified self-loading stretcher - high version



Mercury 7070/4RG PROOF

Stainless steel variable height stretcher certified with 4 rotating wheels



Mercury 7072/4RG PROOF

Certified stainless steel variable height stretcher high version with 4 rotating wheels
