



## 12154

---

"BEAR PODS" Bariatric containment system for self-loading stretchers



---

### FEATURES

**Bear Pods** is designed as an accessory for transporting bariatric patients on self-loading stretchers. The Bariatric Equalizing Abdominal Restraint (BEAR) is made from top quality high-grade materials that are chemical resistant and UV resistant, and may be washed with soap and water. The mesh fabric is constructed of nylon. The purpose of this device is to more safely and reliably transport a bariatric patient in situations involving the use of ambulance stretchers. (Bear Pod fits on all ambulance cots).

This product guarantees you are able to transport a 700 lb. patient on your stretcher, safely and efficiently.

---

### MAIN FEATURES

---


Bariatric



---

## IMAGES

---



---

## ACCESSORY OF

---



**Mercury 7070/4RG PROOF**  
Stainless steel variable height stretcher certified with 4 rotating wheels



**Mercury 7072/4RG PROOF**  
Certified stainless steel variable height stretcher high version with 4 rotating wheels



**Mercury Lite 7080/4RG Proof**  
Ambulance cot in aluminium with variable height and four swivel castors



**Mercury Lite 7082/4RG Proof**  
Certified aluminium variable height stretcher high version with 4 rotating wheels



**MeBer-X 7300**  
"MeBer-X" Multi-level stretcher complete with mattress and belts



**Mercury Cinque 7075/4RG PROOF**  
Certified 5 levels self-loading stretcher with 4 swivel wheels



**Mercury Cinque 7077/4RG PROOF**  
Certified 5 levels self-loading stretcher with 4 swivel wheels - high version



**Mercury Lite Cinque 7095/4RG Proof**  
Certified aluminium 5 levels self-loading stretcher with 4 swivel wheels



**Mercury Lite Cinque 7097/4RG Proof**  
Certified aluminium 5 levels self-loading stretcher with 4 swivel wheels - high version



**Mercury EL 7092 PROOF**  
Certified self-loading stretcher



**Mercury EL 7094 PROOF**  
Certified self-loading stretcher - high version



**7265 Proof**  
Aluminium self loading stretcher "Frog lite Advance" with extractable handles and Trendelenburg