



## 675 Extra Ergolift FX

---

Fixed frame powered Stair Chair suitable for transportation



---

### FEATURES

The range of Me.Ber. powered stair-chairs expands with the new Extra Ergolift FX. Thanks to the great success of the folding model Ergolift, FX borrows to grant the same high performances, but with a fixed frame, especially designed to fulfill the needs of daily transportations services or as a dedicated device in both medical and non medical facilities.

The powered chair is designed for covering the needs of modern all-round service, and it's suitable for transport on vans and ambulances. Ideal for recovering and relocating patients without efforts.

Unique features:

- wireless control panel
- telescopic handle bar
- supplied with thoracic abdominal belt
- long and high gripping tracks
- 4 dedicated points for restraint systems
- 200 mm rear wheels
- high capacity battery
- high comfort

Extra Ergolift loads up to 230 kg and can be used under rain or adverse weather conditions.

Weight 44,8 kg.

28V battery

Autonomy: 18 ramps by 20 steps with a load of 130 kg

Charging time: max 1 hour

MDD 93/42/EEC - Class I Medical Device EN 1789 Certified



---

## MAIN FEATURES


---

Length (cm)	84
Max length (cm)	150
Width (cm)	55
Height (cm)	140
Weight (kg)	44,8
Loading capacity (kg)	230
Cardio	<input checked="" type="checkbox"/>

---

## IMAGES

---



---

## SEE ALSO

---



**Extra Ergolift 671**  
Powered stair chair



**Extra 670/BR**  
Track assisted stair chair with arms



**679**  
Narrow foldable sedan chair

---

---

## ACCESSORIES

---



**12090**  
Battery charger MILWAUKEE M28C  
for rechargeable battery M28 RED  
Extra Ergolift chair



**12092**  
Battery MILWAUKEE M28 RED 5,0  
Ampere for Extra Ergolift chair



**12118**  
Mattress for Extra and Ergolift  
chairs' transformation kit into fixed  
chairs



**988**  
Integrated control fastening  
system for Extra Ergolift FX and  
Proxima Pro chairs (optional)



**566**  
Black bag with 3 belts in 2 pieces  
made of black ribbon with orange  
MeBer logo for Extra chairs



**610**  
Thoracic-abdominal belt with four  
fastening for padded chair



**171J**  
Removable left lever with  
ergonomic handle for Extra chairs.



**172J**  
Removable right lever with  
ergonomic handle for Extra chairs



**670-004**  
Castor with brake diam. 200 mm  
for Extra chair

---

---

## STANDARD ACCESSORIES

---



**12090**  
Battery charger MILWAUKEE M28C  
for rechargeable battery M28 RED  
Extra Ergolift chair



**12092**  
Battery MILWAUKEE M28 RED 5,0  
Ampere for Extra Ergolift chair



**566**  
Black bag with 3 belts in 2 pieces  
made of black ribbon with orange  
MeBer logo for Extra chairs



**12118**  
Mattress for Extra and Ergolift  
chairs' transformation kit into fixed  
chairs



**610**  
Thoracic-abdominal belt with four  
fastening for padded chair

---