



19156-130

"MTX-2R" double space stretcher support with lateral movement and ramp of mm. 130 - EN1789



FEATURES

Double space stretcher support with lateral movement and ramp mm 130 in the underneath space can be stored both a scoop stretcher and a spinal board.
Lateral movement: cm 30

MAIN FEATURES

Length (cm)	205
Width (cm)	62
Weight (kg)	67
Loading capacity (kg)	400
Side sliding	
Plate height	21
Total height	24

IMAGES



SEE ALSO



19156-130 AL

"MTX-2R-AL" double space
stretcher support aluminium alloy
with lateral movement and ramp of
mm. 130



19162

MTX-2-X Double space stainless
steel stretcher support for X-frame
stretchers

ACCESSORIES



19123

Ribaltina 130 mm con attacchi per
supporti e rullo



19115

Ribaltina da 180 mm con attacchi
per supporti e rullo