



## 891/SG

---

Flake basic independent chambers vacuum mattress



---

### FEATURES

Flake basic independent chambers vacuum mattress 200x90.

---

### MAIN FEATURES

Length (cm)	200
Width (cm)	90
Weight (kg)	6,2
Trauma	
Size	Adults
Material	PVC-coated polyester
Loading capacity (kg)	150

---

### SEE ALSO



**"Flake" 891**  
Independent chambers vacuum mattress



**Chesterfield 914**  
Vacuum mattress with independent chambers with central welding