



## 11220 Oxytest

OXYTEST portable pulse oximeter complete with adult sensor, USB cable and AAA battery




### FEATURES

Oxytest is a portable pulse oximeter, easy to use and light for transport. Complete with:

- Oled display, four directions to display arterial oxygen saturation, pulse rate and pulse strength
- battery voltage indicator
- standby or automatic suspension

Thanks to the sensors art. 11224 and art. 11222 (not included in the package), it is possible to use Oxytest also for children and new born.

### MAIN FEATURES

Height (cm)	14,5
Width (cm)	8,3
Depth (cm.)	3,3
Weight (kg)	0,24
Cardio	

### IMAGES



SEE ALSO



**16500**  
**Sterilair Pro**  
Biological air treatment device -  
Wall version



**16502**  
**Sterilair Pro**  
Biological air treatment device -  
Stand version



**Prevention and individual  
protection**



**8377**  
Infrared thermometer



**8506**  
Pulse oximeter Oxy-110



**11200 Vital Test**  
Handheld pulse oximeter  
complete with adult sensor

ACCESSORIES



**11224**  
Pediatric sensor for Oxytest



**11222**  
New born sensor for Oxytest

STANDARD ACCESSORIES



**11226**  
Adult sensor for pulse oximeter  
Oxytest