



## 8432-001

Pedestal for art. 8432



### FEATURES

Cabinet pedestal for art. 8432. Back holes for wall hanging, lower holes for fixing it to the ground and interlocking system for cabinet to be assembled to the pedestal. Made of carbon steel enamelled by epoxy polyester white powder RAL9010 heated at 180°.

### MAIN FEATURES

|             |     |
|-------------|-----|
| Width (cm)  | 25  |
| Height (cm) | 100 |
| Depth (cm.) | 15  |

### SEE ALSO



**8396-001**

Pedestal for art. 8396 – 8398

---

## ACCESSORY OF

---



**8432**

Teca di colore bianco mod.  
standard per defibrillatore

---