



1325/RB

Compact professional rescue backpack Antartide



FEATURES

The “Antartide” model differs from the Artide model in the large front pocket with a crystal window used to transport an AED (any defibrillator model) or other medical material. Made of blue or red nylon and double stitching, large high visibility reflective bands, PVC bottom, reinforced and detachable shoulder straps. The large external pocket enable the rescuer giving a clear view of the contents. Internally it is equipped with Velcro® attachments for fixing three supplied professional bags (blue, yellow and green).

MAIN FEATURES

Height (cm)	44
Width (cm)	36
Depth (cm.)	26
Weight (kg)	3,9

SEE ALSO



1362
“Ambumed” rescue knapsack



1362
“Ambumed” rescue knapsack



1325/RG
Compact professional royal and grey rescue backpack Antartide



1326/RG
Compact professional rescue backpack Artide

STANDARD ACCESSORIES



1362-003/B

Blue professional rectangular inside
bag for knapsack



1362-003/G

Yellow professional rectangular
inside bag for knapsack



1362-003/V

Green professional rectangular
inside bag for knapsack
