



## 8396

---

Defibrillator metal cabinet with thermoregulation and alarm for outer use



---

### FEATURES

Metal cabinet made of p02 carbon steel with 15/10 thickness, watertight welds, gaskets, watertight door and dust with IP44 protection degree and snap closure with tubes adjustment and with the possibility of application of padlock or seal (not supplied). Hinges in durable metal alloy anti tear anti corrosion. Height-adjustable internal shelf with grille for hot air recirculation. Comfortable technopolymer handle for easy opening. Liquid crystal sensor positioned at the front of the temperature inside the case. Alarm signal Volume 100 db.

External dimensions 425 x 480 x 215 mm (W x H x D)  
Internal dimensions 420 mm (W) - 390 mm (H)

ATTENTION: METAL CABINET NOT SUITABLE FOR MARITIME LOCATIONS

---

### SEE ALSO



**8432**  
Teca di colore bianco mod.  
standard per defibrillatore



**8396-001**  
Pedestal for art. 8396 - 8398



**8398**  
Defibrillator metal cabinet with  
thermoregulation and alarm for  
outer use

---

---

## ACCESSORIES

---



**8396-001**

Pedestal for art. 8396 – 8398

---