



8398

Defibrillator metal cabinet with thermoregulation and alarm for outer use



FEATURES

Metal cabinet made of p02 carbon steel with 15/10 thickness, watertight welds, gaskets, watertight door and dust with IP44 protection degree and snap closure with tubes adjustment and with the possibility of application of padlock or seal (not supplied). Hinges in durable metal alloy anti tear anti corrosion. Height-adjustable internal shelf with grille for hot air recirculation. Comfortable technopolymer handle for easy opening. Liquid crystal sensor positioned at the front of the temperature inside the case.

External dimensions 425 x 480 x 215 mm (W x H x D)

Internal dimensions 420 mm (W) - 390 mm (H)

ATTENTION: METAL CABINET NOT SUITABLE FOR MARITIME LOCATIONS

SEE ALSO



8396-001
Pedestal for art. 8396 – 8398



8396
Defibrillator metal cabinet with thermoregulation and alarm for outer use

ACCESSORIES



8396-001

Pedestal for art. 8396 – 8398
