



## 988

---

Integrated control fastening system for Extra Ergolift FX and Proxima Pro chairs (optional)



---

### FEATURES

Quick fixing system to the floor designed to block the padded chair art. 683 and the Extra Ergolift FX chair art. 675 within emergency vehicles and disabled vehicles. Made of stainless steel.

Prepared for integrated release from the backrest.

Device certified according to the European harmonized standards UNI EN 1865- UNI EN 1789.



---

### ACCESSORY OF



**Extra Ergolift FX 675**  
Fixed frame powered Stair Chair  
suitable for transportation



**Proxima Pro 683**  
Fixed padded sedan chair