



19156-130 AL

"MTX-2R-AL" double space stretcher support aluminium alloy with lateral movement and ramp of mm. 130



FEATURES

Same features of art. 19156-130 but with an aluminium structure that reduces the weight of 20 kg.

EN 1789 certified.

Lateral movement: cm 30.

MAIN FEATURES

Width (cm)	62
Length (cm)	205
Weight (kg)	47
Loading capacity (kg)	400
Side sliding	
Plate height	21
Total height	24

SEE ALSO



19156-130

"MTX-2R" double space stretcher support with lateral movement and ramp of mm. 130 - EN1789



19162

MTX-2-X Double space stainless steel stretcher support for X-frame stretchers

ACCESSORIES



19123

Ribaltina 130 mm con attacchi per supporti e rullo



19115

Ribaltina da 180 mm con attacchi per supporti e rullo
