



## 19156-130 AL

"MTX-2R-AL" double space stretcher support aluminium alloy with lateral movement and ramp of mm. 130



### FEATURES

Same features of art. 19156-130 but with an aluminium structure that reduces the weight of 20 kg.

**EN 1789 certified.**

**Lateral movement: cm 30.**

### MAIN FEATURES

Width (cm)	62
Length (cm)	205
Weight (kg)	47
Loading capacity (kg)	400
Side sliding	
Plate height	21
Total height	24

### SEE ALSO



**19156-130**

"MTX-2R" double space stretcher support with lateral movement and ramp of mm. 130 - EN1789



**19162**

MTX-2-X Double space stainless steel stretcher support for X-frame stretchers

---

## ACCESSORIES

---



**19123**

Ribaltina 130 mm con attacchi per supporti e rullo



**19115**

Ribaltina da 180 mm con attacchi per supporti e rullo

---