



## 9114 Chesterfield

---

Vacuum mattress with independent chambers with central welding




---

### FEATURES

CHESTERFIELD vacuum mattress with independent chambers with central welding.  
Usable with X-ray, CT and MRI.

---

### MAIN FEATURES

Trauma	
Length (cm)	210
Width (cm)	100
Weight (kg)	7,8
Loading capacity (kg)	150

---

### IMAGES



---

SEE ALSO

---



**Snake Snow 894/SNOW**  
Independent chambers vacuum  
mattress with anti-slip system for  
snow



**"Flake" 891**  
Independent chambers vacuum  
mattress



**891/SG**  
Flake basic independent chambers  
vacuum mattress



**Snake 894**  
Independent chambers vacuum  
matress



**Naxos 9104**  
"NAXOS" orange and blue vacuum  
mattress without chambers

---