



## 8540

---

Multiparametric monitor TG



---

### FEATURES

The TG monitor is a professional tool with many operational options and memory, fitted with accessories and circuits. The monitor has a dual function: diagnosis and monitoring.

Allows for the individual or simultaneous detection of various vital parameters.

5 parameters


- 5 channel, 2 track ECG
- Systolic/diastolic/average blood pressure
- Heartbeat
- Pulse oximeter
- Breathing control
- Twin probe temperature

Functions

- Multi-parameter monitor with number and track display.
- Operates on battery and mains supply - Power supply: 100 ÷ 240 V ac, 50 ÷ 60 Hz. It is also equipped with a battery to allow continuity of use without external power supply (up to 4 hours with the standard battery)
- Complete with built-in printer
- Network connection for remote management
- Parameter adjustment and setting
- 10" monitor

---

### MAIN FEATURES

Cardio	
Height (cm)	24,6
Width (cm)	26,1
Depth (cm.)	14,6
Weight (kg)	2,8

---

## ACCESSORIES

---



**8540-001**

Reusable Spo2 sensor for adults for multi-parameter monitor



**8540-002**

27-35 cm adult cuff



**8540-005**

Rechargeable lithium battery



**8540-006**

Thermal graph paper



**8540-003**

NIBP bracelet, E7, pediatric, 16-21.5 cm



**8540-004**

NIBP cuff, E5, neonatal, 10-15 cm

---

## STANDARD ACCESSORIES

---



**8540-001**

Reusable Spo2 sensor for adults for multi-parameter monitor



**8540-002**

27-35 cm adult cuff