



8540

Multiparametric monitor TG



FEATURES

The TG monitor is a professional tool with many operational options and memory, fitted with accessories and circuits. The monitor has a dual function: diagnosis and monitoring.

Allows for the individual or simultaneous detection of various vital parameters.


5 parameters

- 5 channel, 2 track ECG
- Systolic/diastolic/average blood pressure
- Heartbeat
- Pulse oximeter
- Breathing control
- Twin probe temperature

Functions

- Multi-parameter monitor with number and track display.
- Operates on battery and mains supply - Power supply: 100 ÷ 240 V ac, 50 ÷ 60 Hz. It is also equipped with a battery to allow continuity of use without external power supply (up to 4 hours with the standard battery)
- Complete with built-in printer
- Network connection for remote management
- Parameter adjustment and setting
- 10" monitor

MAIN FEATURES

Cardio	
Height (cm)	24,6
Width (cm)	26,1
Depth (cm.)	14,6
Weight (kg)	2,8

ACCESSORIES



8540-001

Reusable Spo2 sensor for adults for multi-parameter monitor



8540-002

27-35 cm adult cuff



8540-005

Rechargeable lithium battery



8540-006

Thermal graph paper



8540-003

NIBP bracelet, E7, pediatric, 16-21.5 cm



8540-004

NIBP cuff, E5, neonatal, 10-15 cm

STANDARD ACCESSORIES



8540-001

Reusable Spo2 sensor for adults for multi-parameter monitor



8540-002

27-35 cm adult cuff
