



11023

10g certified wall bracket with integrated charging system for CORMOS and CORMOS PLUS



FEATURES

10g certified wall bracket with integrated charging system for CORMOS and CORMOS PLUS.

Dimensions: about 240 x 255 x 5 mm

MAIN FEATURES

Weight (kg)

0,8

ACCESSORY OF



Cormos 11010/FA
Portable aspirator for secretions
with autoclavable canister



Cormos 11010/LN
Portable aspirator for secretions
with disposable liner



Cormos Plus 11022/FA
White complete medical suction
unit with autoclavable canister



Cormos Plus 11022/LN
White complete medical suction
unit with canister for disposable
unit