



707

Red 6-handle fireproof cover with Me.Ber tape for self-loading stretcher



FEATURES

Slash is MeBer's innovative mattress protector transfer sheet. Excellent ally for moving the patient from bed to stretcher and vice versa. Equipped with 6 handles and made of fireproof, tearproof and waterproof nylon. Placed on self-loading stretchers, it easily lends itself to immediate sanitation between one patient and another. The reduced sizes are specially designed to facilitate operations during the handling phases. Compatible with all stretcher models.

SLASH IS NOT A CARRYING SHEET.

MAIN FEATURES

Length (cm)	188
Width (cm)	50
Weight (kg)	1
Loading capacity (kg)	160

IMAGES



ACCESSORY OF



Mercury Lite 7082/4RG Proof
Certified aluminium variable height
stretcher high version with 4
rotating wheels



Mercury Cinque 7075/4RG PROOF
Certified 5 levels self-loading
stretcher with 4 swivel wheels



Mercury Cinque 7077/4RG PROOF
Certified 5 levels self-loading
stretcher with 4 swivel wheels -
high version



**Mercury Lite Cinque 7095/4RG
Proof**
Certified aluminium 5 levels self-
loading stretcher with 4 swivel
wheels



**Mercury Lite Cinque 7097/4RG
Proof**
Certified aluminium 5 levels self-
loading stretcher with 4 swivel
wheels - high version



Mercury Lite 7080/4RG Proof
Ambulance cot in aluminium with
variable height and four swivel
castors



Mercury 7070 PROOF
Certified stainless steel variable
height stretcher



Mercury 7070/4RG PROOF
Stainless steel variable height
stretcher certified with 4 rotating
wheels



Mercury 7072 PROOF
Certified stainless steel variable
height stretcher high version



Mercury 7072/4RG PROOF
Certified stainless steel variable
height stretcher high version with 4
rotating wheels



Mercury EL 7092 PROOF
Certified self-loading stretcher



Mercury EL 7094 PROOF
Certified self-loading stretcher -
high version



7265 Proof
Aluminium self loading stretcher
"Frog lite Advance" with
extractable handles and
Trendelenburg
