



668 Kuma

Black powered stair chair



FEATURES

Kuma is the entry level model of the MeBer powered chairs. Characterized by aluminum alloy frame, lithium-ion battery, front levers, practical armrests, footrest, practical extendable handle (with attached small console for controlling the tracks) and 4 wheels (2 fixed Ø200 mm with brake and 2 swiveling Ø125 mm). The electric motor drives the tracks allowing the ascent and descent of the stairs without effort for the operator. Kuma guarantees a payload of 181 kg and offers a comfortable seat thanks to the standard cushion and headrest. Practical and easy to handle, the Kuma motorized chair is very simple to use. Black colour.

Available with the following standard accessories:

- Seat cushion
- Headrest
- Battery
- Charger
- Belts

MAIN FEATURES

Length (cm)	79
Width (cm)	52
Min height (cm)	108
Max height (cm)	164
Epaisseur pliée (cm)	24,2
Weight (kg)	34
Loading capacity (kg)	181

ACCESSORIES



977
Fixing system for Kuma stair chair



12182
Battery charger for KUMA
motorized chair



12180
Battery for Kuma motorized chair

STANDARD ACCESSORIES



12180
Battery for Kuma motorized chair



12182
Battery charger for KUMA
motorized chair