



## 668 Kuma

---

Black powered stair chair



---

### FEATURES

Kuma is the MeBer powered chair. Characterized by aluminum alloy frame, lithium-ion battery, front levers, practical armrests, footrest, practical extendable handle (with attached small console for controlling the tracks) and 4 wheels (2 fixed Ø200 mm with brake and 2 swiveling Ø125 mm). The electric motor drives the tracks allowing the ascent and descent of the stairs without effort for the operator. Kuma guarantees a payload of 181 kg and offers a comfortable seat thanks to the standard cushion and headrest. Practical and easy to handle, the Kuma motorized chair is very simple to use. Black colour.

Available with the following standard accessories:

- Seat cushion
- Headrest
- Battery
- Charger
- Belts

---

### MAIN FEATURES

Length (cm)	79
Width (cm)	52
Min height (cm)	108
Max height (cm)	164
Epaisseur pliée (cm)	24,2
Weight (kg)	34
Loading capacity (kg)	181

---

## ACCESSORIES

---



**977**  
Fixing system for Kuma stair chair



**12182**  
Battery charger for KUMA  
motorized chair



**12180**  
Battery for Kuma motorized chair

---

## STANDARD ACCESSORIES

---



**12180**  
Battery for Kuma motorized chair



**12182**  
Battery charger for KUMA  
motorized chair