



# 7094 PROOF Mercury EL

Certified self-loading stretcher - high version



### **FEATURES**

Mercury EL is the newest entry in Me.Ber's range of ambulance cots. Sharing many features from other Mercury models, the EL is a naturally highly reliable and versatile cot. This model represents the entry level for the Mercury line.

Unlike most Mercury models, the EL does not feature a foldable front trolley nor 4 swivel castors and it's not multilevel. However, this version serves as a **perfect solution** when a high standard of quality must be met while keeping things **easy and accessible**.

#### Safe-Frame Technology:

The frame of an ambulance cot must withstand intense mechanical stress: this is why Mercury features a single-piece frame, in line with Me.Ber.'s Same Frame building philosophy.

Thanks to this expedient, the structure of the cot is extremely safe and robust through time.



Device certified according to the European harmonized safety standards UNI EN 1865 - UNI EN 1789.

\*All Pictures shown are for illustration purpose only.

# MAIN FEATURES

Type	
Size	Adults
Usage	Transport
10G Certification	<b>⊌</b>
Wheel type	2 Rotating Wheels + 2 Big Fixed Wheels
Lying surface	With Trendelenburg
Max height (cm)	99
Max loading height (cm)	71
Min loading height (cm)	69
Folded height (cm)	31
Width (cm)	57
Length (cm)	197
Weight (kg)	46,6
Loading capacity (kg)	250
Material	Painted steel ABS HDPE (High-density polyethylene)

### **SEE ALSO**



Mercury EL 7092 PROOF Certified self-loading stretcher



Mercury 7070/4RG PROOF Stainless steel variable height stretcher certified with 4 rotating wheels



Mercury 7072/4RG PROOF Certified stainless steel variable height stretcher high version with 4 rotating wheels



Mercury Lite 7080/4RG Proof Ambulance cot in aluminium with variable height and four swivel castors



Mercury Lite 7082/4RG Proof Certified aluminium variable height stretcher high version with 4 rotating wheels



Mercury Cinque 7075/4RG PROOF Certified 5 levels self-loading stretcher with 4 swivel wheels



Mercury Cinque 7077/4RG PROOF Certified 5 levels self-loading stretcher with 4 swivel wheels high version



Mercury Lite Cinque 7095/4RG Proof Certified aluminium 5 levels selfloading stretcher with 4 swivel wheels



Mercury Lite Cinque 7097/4RG Proof Certified aluminium 5 levels selfloading stretcher with 4 swivel wheels - high version

### **ACCESSORIES**



12025
Black multi usage object holder for automatic loading stretcher



12154
"BEAR PODS" Bariatric
containment system for selfloading stretchers



**12018N**Telescopic foldable phleboclysis stand for self loading stretchers



939
Instrument support for stretchers



12080
"MPB" wide bariatric platform for self-loading stretchers Frog
Mercury and Gecko complete of belts with reel and fixing for stretcher



**12080-001**Monting kit for bariatric platform on self loading stretcher



**12084-001** Cover for bariatric wings set



**12070/BG**"COVER OX 2-3 BG"- professional cylinder holder for 2 and 3 lt oxygen bottles



12070/BR COVER OX2-3 BR professional cylinder holder with belt and orange MeBer logo for 2 and 3 lt oxygen bottles



**12072/BG**"COVER OX 5 BG"- professional cylinder holder for 5 lt oxygen bottles



**12072/BR**"COVER OX 5 BR"- professional cylinder holder for 5 lt oxygen bottles



955 PROOF EN 1789 fixing system for Proof rollin stretchers



Mattress for AIR-COM stretchers -Anti-bedsore mattress with alternate segments



**7013/A** 3 pcs. Skay fire-proof electrowelded mattress



**7013/N** 3 pcs. Skay fire-proof electrowelded mattress



938 Oximeb oxygen cylinder holder



**608/A-MEB**2 pieces orange belt with metallic hook and MeBer band (L. 155 cm)



**608/N-MEB**2 pieces black belt with metalic hook and MeBer band (L. 155 cm)



**612/AN-2/MEB**Adjustable thoracic abdominal belt



**612/N-2/MEB** Adjustable thoracic abdominal belt



Bariatric wings set for Mercury self-loading stretchers complete with mattresses, belts with reel and assembly kit



**707**Red 6-handle fireproof cover with
Me.Ber tape for self-loading
stretcher

# STANDARD ACCESSORIES



**608/N-MEB**2 pieces black belt with metalic hook and MeBer band (L. 155 cm)



612/N-2/MEB

Adjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self loading stretchers and sedan chairs



7013/N 3 pcs. Skay fire-proof electrowelded mattress