



## 19170

Double-space ambulance space cot support in aluminium, with reduced height, lateral movement, HERCULES system, 955 and 180 mm ramp with roller




### FEATURES

The HERCULES is the new ambulance cot support with motor-powered loading and unloading. The integrated position detectors and small pushbutton ensure a safe and easy use of this support. The loading takes place by means of a quick-release buckle system with belt, actioned by an integrated coaxial motor. No outer volume. Equipped with a double-safety system, to ensure total safety during loading and unloading). Can be installed on MTX-2R and MTX-2R-AL. 180 mm ramp and 955 fixing system. The underneath space may house both a scoop stretcher and a spinal board, despite the reduced height (-40 cm) compared to the MTX-2.

**Certified EN1789.**

**30 cm lateral movement.**

### MAIN FEATURES

Bariatric	
Length (cm)	216,5
Width (cm)	57,8
Weight (kg)	65
Loading capacity (kg)	400
Plate height	21
Total height	34

---

## IMAGES

---



---

## ACCESSORIES

---



**19123**

Ribaltina 130 mm con attacchi per supporti e rullo



**19115**

Ribaltina da 180 mm con attacchi per supporti e rullo

---