



7265 Proof

Aluminium self loading stretcher "Frog lite Advance" with extractable handles and Trendelenburg



FEATURES

Entirely made of painted aluminum alloy, very light, it is also equipped with trendelenburg. 4 polyurethane wheels (2 fixed with brakes and 2 self-centering swivel Ø 200 mm) offer high comfort on uneven floors and irregular surfaces. The material with which it is built makes it easy to sanitize.

Its main features are:

- Front and rear extractable levers to facilitate rescuers' manoeuvres
- Rigid platform in painted aluminum alloy
- 0°-80° reclining backrest lockable in any position
- Folding side rails
- Vertical legs protected by bumpers
- Swivel wheels equipped with self-centering system
- Folding front landing gear



Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789

MAIN FEATURES

Length (cm)	196
Length with folded trolley (cm)	167
Width (cm)	56
Max height (cm)	91,5
Folded height (cm)	31
Min loading height (cm)	60
Max loading height (cm)	66
Weight (kg)	38,3
Loading capacity (kg)	250

IMAGES



SEE ALSO



Frog Lite 7260 PROOF
"FROG LITE" self loading aluminium
stretcher with trendelenburg
certified



Frog 7240 PROOF
"FROG PLUS" self loading stretcher
with trendelenburg certified



Frog 7210 PROOF
"FROG" yellow self loading stretcher
with trendelenburg certified

ACCESSORIES



939
Instrument support for stretchers



936
Phleboclysis stand for automatic loading stretcher Winner and Frog



955 PROOF
EN 1789 fixing system for Proof roll-in stretchers



707
Red 6-handle fireproof cover with Me.Ber tape for self-loading stretcher



12080
"MPB" wide bariatric platform for self-loading stretchers Frog Mercury and Gecko complete of belts with reel and fixing for stretcher



12084
Bariatric wings set for self-loading stretchers complete with mattresses and belts with reel



938
Oximeb oxygen cylinder holder



12054
Oximeb Maxi - Oxygen cylinder holder



12128
Adjustable system for headrest of Meber-X and Mercury stretchers



12106
Mattress for AIR-COM stretchers - Anti-bedsores mattress with alternate segments



12154
"BEAR PODS" Bariatric containment system for self-loading stretchers



12070/BG
"COVER OX 2-3 BG" - professional cylinder holder for 2 and 3 lt oxygen bottles



12070/BR
COVER OX 2-3 BR professional cylinder holder with belt and orange MeBer logo for 2 and 3 lt oxygen bottles



12072/BG
"COVER OX 5 BG" - professional cylinder holder for 5 lt oxygen bottles



12072/BR
"COVER OX 5 BR" - professional cylinder holder for 5 lt oxygen bottles



7013/N
3 pcs. Skay fire-proof electrowelded mattress



608/N-MEB
2 pieces black belt with metallic



hook and MeBer band (L. 155 cm)

608/A-MEB

2 pieces orange belt with metallic hook and MeBer band (L. 155 cm)



612/AN-2/MEB

Adjustable thoracic abdominal belt orange/black with MeBer band and 4 points fastening for self loading stretchers and sedan chairs



612/N-2/MEB

Adjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self loading stretchers and sedan chairs

STANDARD ACCESSORIES



608/N-MEB

2 pieces black belt with metallic hook and MeBer band (L. 155 cm)



612/N-2/MEB

Adjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self loading stretchers and sedan chairs



7013/N

3 pcs. Skay fire-proof electrowelded mattress