



## 7265 Proof

---

Aluminium self loading stretcher "Frog lite Advance" with extractable handles and Trendelenburg



---

### FEATURES

Entirely made of painted aluminum alloy, very light, it is also equipped with trendelenburg. 4 polyurethane wheels (2 fixed with brakes and 2 self-centering swivel Ø 200 mm) offer high comfort on uneven floors and irregular surfaces. The material with which it is built makes it easy to sanitize.

Its main features are:

- Front and rear extractable levers to facilitate rescuers' manoeuvres
- Rigid platform in painted aluminum alloy
- 0°-80° reclining backrest lockable in any position
- Folding side rails
- Vertical legs protected by bumpers
- Swivel wheels equipped with self-centering system
- Folding front landing gear



Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789

---

## MAIN FEATURES

---

Length (cm)	196
Length with folded trolley (cm)	167
Width (cm)	56
Max height (cm)	91,5
Folded height (cm)	31
Min loading height (cm)	60
Max loading height (cm)	66
Weight (kg)	38,3
Loading capacity (kg)	250

---

## IMAGES

---



**Frog Lite 7260 PROOF**  
"FROG LITE" self loading aluminium  
stretcher with trendelenburg  
certified



**Frog 7240 PROOF**  
"FROG PLUS" self loading stretcher  
with trendelenburg certified



**Frog 7210 PROOF**  
"FROG" yellow self loading stretcher  
with trendelenburg certified

---

## ACCESSORIES



**939**  
Instrument support for stretchers



**936**  
Phleboclysis stand for automatic loading stretcher Winner and Frog



**955 PROOF**  
EN 1789 fixing system for Proof roll-in stretchers



**707**  
Red 6-handle fireproof cover with Me.Ber tape for self-loading stretcher



**12080**  
"MPB" wide bariatric platform for self-loading stretchers Frog Mercury and Gecko complete of belts with reel and fixing for stretcher



**12084**  
Bariatric wings set for self-loading stretchers complete with mattresses and belts with reel



**938**  
Oximeb oxygen cylinder holder



**12054**  
Oximeb Maxi - Oxygen cylinder holder



**12128**  
Adjustable system for headrest of Meber-X and Mercury stretchers



**12106**  
Mattress for AIR-COM stretchers - Anti-bedsores mattress with alternate segments



**12154**  
"BEAR PODS" Bariatric containment system for self-loading stretchers



**12070/BG**  
"COVER OX 2-3 BG"- professional cylinder holder for 2 and 3 lt oxygen bottles



**12070/BR**  
COVER OX2-3 BR professional cylinder holder with belt and orange MeBer logo for 2 and 3 lt oxygen bottles



**12072/BG**  
"COVER OX 5 BG"- professional cylinder holder for 5 lt oxygen bottles



**12072/BR**  
"COVER OX 5 BR"- professional cylinder holder for 5 lt oxygen bottles



**7013/N**  
3 pcs. Skay fire-proof electrowelded mattress



**608/N-MEB**  
2 pieces black belt with metallic



hook and MeBer band (L. 155 cm)

**608/A-MEB**

2 pieces orange belt with metallic hook and MeBer band (L. 155 cm)



**612/AN-2/MEB**

Adjustable thoracic abdominal belt orange/black with MeBer band and 4 points fastening for self loading stretchers and sedan chairs



**612/N-2/MEB**

Adjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self loading stretchers and sedan chairs

**STANDARD ACCESSORIES**



**608/N-MEB**

2 pieces black belt with metallic hook and MeBer band (L. 155 cm)



**612/N-2/MEB**

Adjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self loading stretchers and sedan chairs



**7013/N**

3 pcs. Skay fire-proof electrowelded mattress