




12200

Bariatric wings set for Mercury self-loading stretchers complete with mattresses, belts with reel and assembly kit



FEATURES

The bariatric wings are a practical accessory that allows you to convert the self-loading stretcher into a practical bariatric recovery tool. Thanks to the small footprint (and a weight of only 14 kg) they can be easily transported by ambulance to be used when needed. Their use allows to face in safety, and in the most adequate way possible, every recovery scenario in the presence of bariatric patients. The kit includes two wings that extend the support surface of the stretcher by over 40 cm; assembly and disassembly on the self-loading stretcher are extremely simple and intuitive and do not require the use of tools. The bariatric wings are equipped with mattresses (art. 12170 for right bariatric wing and art. 12172 for left bariatric wing), belts with reel and assembly kit art. 12080-001. They can be purchased in kits including both the right and the left (art. 12200) or separately (art. 12202 Right bariatric wing and 12204 left bariatric wing).


For Frog and MBS-01 ambulance cots use art. [12084](#).

Dimensions:
1365x205x85 mm (a wing)

MAIN FEATURES

Length (cm)	136,5
Width (cm)	20,5

IMAGES



SEE ALSO



12084

Bariatric wings set for self-loading stretchers complete with mattresses and belts with reel

ACCESSORIES



12170

Black mattress for right bariatric wing for self-loading stretchers



12172

Black mattress for left bariatric wing for self-loading stretchers

STANDARD ACCESSORIES



12080-001

Monting kit for bariatric platform on self loading stretcher

ACCESSORY OF



Mercury Lite 7080/4RG Proof
Ambulance cot in aluminium with
variable height and four swivel
castors



Mercury 7070/4RG PROOF
Stainless steel variable height
stretcher certified with 4 rotating
wheels



Mercury 7072/4RG PROOF
Certified stainless steel variable
height stretcher high version with 4
rotating wheels



Mercury Lite 7082/4RG Proof
Certified aluminium variable
height stretcher high version with
4 rotating wheels



Mercury Cinque 7075/4RG PROOF
Certified 5 levels self-loading
stretcher with 4 swivel wheels



Mercury Cinque 7077/4RG PROOF
Certified 5 levels self-loading
stretcher with 4 swivel wheels -
high version



**Mercury Lite Cinque 7095/4RG
Proof**
Certified aluminium 5 levels self-
loading stretcher with 4 swivel
wheels



**Mercury Lite Cinque 7097/4RG
Proof**
Certified aluminium 5 levels self-
loading stretcher with 4 swivel
wheels - high version



Mercury EL 7092 PROOF
Certified self-loading stretcher



Mercury EL 7094 PROOF
Certified self-loading stretcher -
high version