



# 629-007

Slim black plate for art. 629 head immobilizer



#### **FEATURES**

The SLIM Black Base is a low-profile PVC board designed to facilitate spine board storage in the ambulance holder, even with the head immobilizer installed.

Thanks to its slim design, it effectively solves the space issues encountered with traditional bases.

Made of **high-strength PVC**, the base is **easy to clean and disinfect**, ensuring long-lasting durability. It comes **complete with straps** for secure and stable attachment of the head immobilizer to the spine board.

Low-profile design for easier storage Material: hygienic, easy-to-clean PVC Straps included

## MAIN FEATURES

Width (cm)	38,5
Depth (cm.)	26
Thickness (cm)	1,2

#### **IMAGES**



## **ACCESSORY OF**



**Jelly 629** Red standard head immobilizer for spinal board