



## 629-007

---

Slim black plate for art. 629 head immobilizer



---

### FEATURES

The **SLIM Black Base** is a **low-profile PVC board** designed to facilitate **spine board storage in the ambulance holder**, even with the head immobilizer installed.

Thanks to its **slim design**, it effectively solves the space issues encountered with traditional bases.

Made of **high-strength PVC**, the base is **easy to clean and disinfect**, ensuring long-lasting durability. It comes **complete with straps** for secure and stable attachment of the head immobilizer to the spine board.

Low-profile design for easier storage

Material: hygienic, easy-to-clean PVC

Straps included

---

### MAIN FEATURES

Width (cm)	38,5
Depth (cm.)	26
Thickness (cm)	1,2

---

## IMAGES

---



---

## SEE ALSO

---



**629/SLIM**  
Red head immobilizer with black  
slim base for spinal board

---

---

## ACCESSORY OF

---



**Jelly 629**  
Red standard head immobilizer for  
spinal board

---