



8548

Multiparameter monitor PC-300 SpO2, NIBP



FEATURES

PC-300 handheld Spot-check monitor measures SpO2, Pulse Rate, NIBP and TEMP. Widely used in hospitals and clinics. Equipped with Bluetooth and USB interface, allows to transfer data to PC, tablets and smartphone (Android and iOS).

- up to 6000 groups NIBP data storage, up to 100 different patients (each user can store 999 measurement records)
- quick ear probe TEMP measurement result in 5 seconds
- data transfer to PC/Smartphone by USB or Bluetooth
- optional NIBP pediatric cuff and SpO2 pediatric probe
- optional ECG and Glucometer
- built-in rechargeable lithium battery
- communication protocols available on request

Power Supply

Battery: 3.7 V built-in Li-ion battery

Operating voltage: 3.6-4.2 VDC

Operating current: ≈800 mA

Charger: 100 - 240 VAC, 50/60 Hz, 5 VDC or USB

Supplied with:

SpO2 adult probe

NIBP adult cuff


Ear thermometer (probe covers not included)

Lithium battery with charger

Carrying bag

User manual

MAIN FEATURES

Cardio	
Length (cm)	16,5
Thickness (cm)	6,8
Width (cm)	9,6

IMAGES



SEE ALSO



B1 CLASSIC 8534
Automatic blood pressure with
adult cuff

ACCESSORIES

8548-001
Sensore SpO2 pediatrico per
monitor multiparametrico PC-300
art. 8548



8548-002
Sensore SpO2 adulti per monitor
multiparametrico PC-300 art. 8548

8548-003
Bracciale pediatrico 10-19 cm per
monitor multiparametrico PC-300



8548-004
Bracciale adulto per monitor
multiparametrico PC-300

8548-005
Bracciale bariatrico per monitor
multiparametrico PC-300