




1550 Folding stretcher

Military folding stretcher with green canvas



FEATURES

Width-collapsible stretcher built in tubular aluminium and green fireproof military canvas. It features telescopic handles and strap belts for patient immobilization. Ideal for military applications.

MAIN FEATURES

Size	Adults
Usage	Transport
Foldable	
Lying surface	Without Trendelenburg
Height (cm)	17
Width (cm)	58
Length (cm)	198
Max length (cm)	231
Epaisseur pliée (cm)	8,5
Weight (kg)	9
Loading capacity (kg)	170
Material	Anodized aluminium Nylon

IMAGES



SEE ALSO



Folding stretcher 660
Collapsible stretcher foldable in 4
pieces, telescopic handles

STANDARD ACCESSORIES



605/Nmeb
50 mm black belt in one piece with
plastic buckle and ME.BER. ribbon
