



16100 Toboga



Stretcher with straps and foot rest



FEATURES

Consisting of an ABS modified polymer in the outer shell, supported by a strong aluminium frame. 12 hand grips ensure a very easy handling of the stretcher. Chemical-, UV-, rust- and corrosion resistant. Robust brass grommets, positioned on the sides, enable lifting the stretcher by means of art. 16100-005. One layer of EVA foam ensures high comfort to the patient. Supplied with 4 straps and adjustable foot rest.

MAIN FEATURES

Size	Adults
Usage	Rescue Immobilization Transport
Height (cm)	17
Width (cm)	61
Length (cm)	216
Weight (kg)	15
Loading capacity (kg)	300
Trauma	
Alpine	
Material	ABS Aluminium alloy

SEE ALSO



Toboga 16102
Detachable with straps and foot
rest

ACCESSORIES



16100-002
Flotation system for toboga
stretcher



16100-005
Certified lifting sling for Toboga
stretchers
