



16102 Toboga

Detachable with straps and foot rest



FEATURES

Consisting of an ABS modified polymer in the outer shell, supported by a strong aluminium frame. 12 hand grips ensure a very easy handling of the stretcher. Chemical-, UV-, rust- and corrosion resistant. Robust brass grommets, positioned on the sides, enable lifting the stretcher by means of art. 16100-005. One layer of EVA foam ensures high comfort to the patient. Supplied with 4 straps and adjustable foot rest.

This model is made up of 2 halves that enable easy storage and carrying. The two halves are kept together by strong aluminium rods locked together with steel pivots.

MAIN FEATURES

| Height (cm) | 17 |
|-----------------------|---------------------------------------|
| Width (cm) | 61 |
| Length (cm) | 216 |
| Weight (kg) | 15 |
| Loading capacity (kg) | 300 |
| Size | Adults |
| Usage | Rescue Immobilization Transport |
| Trauma | ₩ |
| Alpine | ₩ |
| Material | ABS Aluminium alloy |

IMAGES



SEE ALSO



ACCESSORIES







16100-005 Certified lifting sling for Toboga stretchers