



## 16102 Toboga

Detachable with straps and foot rest





### FEATURES

Consisting of an ABS modified polymer in the outer shell, supported by a strong aluminium frame. 12 hand grips ensure a very easy handling of the stretcher. Chemical-, UV-, rust- and corrosion resistant. Robust brass grommets, positioned on the sides, enable lifting the stretcher by means of art. 16100-005. One layer of EVA foam ensures high comfort to the patient. Supplied with 4 straps and adjustable foot rest.

This model is made up of 2 halves that enable easy storage and carrying. The two halves are kept together by strong aluminium rods locked together with steel pivots.

### MAIN FEATURES

Height (cm)	17
Width (cm)	61
Length (cm)	216
Weight (kg)	15
Loading capacity (kg)	300
Size	Adults
Usage	Rescue Immobilization Transport
Trauma	
Alpine	
Material	ABS Aluminium alloy

---

IMAGES

---



---

SEE ALSO

---



**Toboga 16100**  
Stretcher with straps and foot rest

---

ACCESSORIES

---



**16102-001**  
Carry bag for detachable Toboga  
stretcher



**16100-002**  
Flotation system for toboga  
stretcher



**16100-005**  
Certified lifting sling for Toboga  
stretchers