



16102 Toboga

Detachable with straps and foot rest



FEATURES

Consisting of an ABS modified polymer in the outer shell, supported by a strong aluminium frame. 12 hand grips ensure a very easy handling of the stretcher. Chemical-, UV-, rust- and corrosion resistant. Robust brass grommets, positioned on the sides, enable lifting the stretcher by means of art. 16100-005. One layer of EVA foam ensures high comfort to the patient. Supplied with 4 straps and adjustable foot rest.

This model is made up of 2 halves that enable easy storage and carrying. The two halves are kept together by strong aluminium rods locked together with steel pivots.

MAIN FEATURES

Height (cm)	17
Width (cm)	61
Length (cm)	216
Weight (kg)	15
Loading capacity (kg)	300
Size	Adults
Usage	Rescue Immobilization Transport
Trauma	
Alpine	
Material	ABS Aluminium alloy

IMAGES



SEE ALSO



Toboga 16100
Stretcher with straps and foot rest

ACCESSORIES



16102-001
Carry bag for detachable Toboga
stretcher



16100-002
Flotation system for toboga
stretcher



16100-005
Certified lifting sling for Toboga
stretchers
