



## 1052

Oxygen therapy mask for adults at variable concentration - venturi system



### FEATURES

Soft oxygen therapy mask at variable concentration Venturi System complete with coloured connections kit to change concentration of given oxygen. Anti-flattening connection tube.

### MAIN FEATURES

Size

Adults

### SEE ALSO



**1055**

Oxygen therapy mask for adults with bag



**1050**

Oxygen therapy mask for adults



**1045**

Oxygen therapy mask for children with bag



**1040**

Oxygen therapy mask for children



**1053**

Adult mask for inhalation therapy

---

## ACCESSORIES

---



**1057**

Goggles for oxygen therapy with rear sliding fastener



**1056**

Extension tube for oxygen f/f 2 metres long

---