



1740

Rotating wall aneroid sphygmomanometer



FEATURES

Round sphygmomanometer. Easy to read, it is ideal for ambulances and ambulatories. 65° rotating in both directions, roomy cuff basket for storage of accessories, coiled tube extendible to a length of 3 m, latex bulb and chrome-plated air release valve, highly legible dial-plate diam. 147 mm. Specially tempered copper-beryllium diaphragm. Calibrated nylon Velcro[®] cuff with 2-tube latex bladder.

IMAGES



ACCESSORIES



1740-001

Pediatric cuff for rotating wall
aneroid sphygmomanometer Art.
1740.
