



1880-002

Mattress pillow to be used with Icarus stretcher



FEATURES

Mattress pillow to be used with Icarus stretcher. Made of inner fireproof padding, waterproof, anti-tear Skay cover. Double seam.

Dimensions: cm 55x35x5

*All Pictures shown are for illustration purpose only.

MAIN FEATURES

Width (cm)	55
Length (cm)	35
Thickness (cm)	5
Material	Skay

IMAGES



SEE ALSO



1880-001
Mattress for Icarus

ACCESSORIES



Icarus 13010
Multi level stretcher with side rails

ACCESSORY OF



1880-001
Mattress for Icarus
